




PrismView™ Blades



Laryngoscopes With *FIBER OPTIC* Illumination.



End View

The NEW PrismView™ Fiber Optic blades are available in your choice of an American or English Profile. PrismView™ blades may be used: a) when the view of the vocal chords is impaired, especially the most anterior larynx. b) during nasal intubation and c) for post operative examination of the larynx. Both contain an optically polished prism which provides a refraction of 30° without image inversion. To prevent condensation, warm the prism before use by immersing the complete unit in warm water.

PLEASE NOTE: The Sun•Med GreenLine® Fiber Optic System is COMPATIBLE with the following companies “GREEN” Systems: Welch Allyn, Rüsch, Propper, Heine, AMS, and Amtec



ENGLISH	AMERICAN	SIZE	DESCRIPTION	LENGTH	DEPTH
5-5228-02	5-5229-02	2	Child	100mm	14mm
5-5228-03	5-5229-03	3	Medium Adult	130mm	24mm
5-5228-35	5-5229-35	3½	Med Adult Extended	144mm	24mm
5-5228-04	5-5229-04	4	Large Adult	155mm	24mm